



TePe Gentle Care™

Handle Polypropylene (ISCC certified).

Filaments Polyamide (PA).
TePe Gentle Care has 5400 super soft filaments.

Pack *Consumer pack; 1 pc/bag*
Retail carton; 25 pcs/carton
The bag is made of Polypropylene (PP).

Consumer pack; 1 pc/blister pack
Retail carton; 14 pcs/carton
The blister is made of polyethylene terephthalate, (PET), partly recycled.
The backing card is made of FSC certified cardboard.
The tray in the retail box is made of corrugated cardboard.

Retail and transport cartons for all packaging above; Chlorine-free cardboard respectively corrugated cardboard. All from recycled material, except for the retail cartons which are from partly recycled material and partly from virgin FSC certified material.

Recycling Toothbrush in household waste.

Bag in container for plastic recycling.

The blister pack should be separated and recycled; the blister in container for plastic recycling and the backing card in container for recycling of cardboard.

The cartons in container for recycling of cardboard/corrugated cardboard.

Product features

- Super soft filaments for gentle cleaning.
- Second brush stage after oral surgery.
- For patients with dry mouth, sensitive root surfaces and sensitive oral tissues.
- Neck can be angled for easier reach.

Examples of areas of use TePe Gentle Care is a super soft toothbrush recommended for brushing after oral surgery and for patients with sore or very delicate oral tissue, very sensitive root surfaces, extremely dry mouth, and after radiotherapy. Additional areas of use can occur upon recommendation from a dental professional. Consult your dental professional regarding the use of TePe Gentle Care.

Developed in collaboration with dental experts Prof. and Chairman Hiroshi Okamoto
Faculty of Dentistry, Ohu University,
Tokyo, Japan.